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ENGLISH

1. Remove the covers from the back wheel. Remove the rubber plugs from the outside of the brake wheel with a screwdriver, and pull out the hexagonal axle by using an 8 mm bolt. Remove the brake shaft and brake stick (A and instruction 105869).
2. Drill an 8.5 mm hole through the top plate of the underframe. Position the hole in the centre, approximately 7 cm from the outer edge on the right-hand side (B).
3. Mount the fitting for the gas cartridge under the frame with a 6 mm bolt in the two holes with the clasp on the top side (C and D).
4. Mount the brake shaft and hexagonal axle through the wheel and brake shaft. Lock the axle by mounting the two rubber plugs on the outside of the wheels once again.
5. Mount the actuator with bolts and clasps (kit). Assemble the plug from spiral cable and push it in under the fitting (C).
6. Mount the rubber feed plug in the hole on the frame, and insert the spiral wire with mox plug through the hole (B and C).
7. Remove the battery covers and connect the motor control to the batteries. Clip the existing plug off the wires, and connect and clamp the two wires together in the cable shoes. NB: Connect the red wire with red to the plus pole, and connect the black wire with black to the minus pole. Re-attach the battery covers (F).
8. Connect the Mini-Fit plug to the spiral cable, and connect this cable to the other Mini-Fit plug from the motor control (G and H).
9. Pull up the spiral cable from the motor control to the armrest cushion, and set the jack plug in the bushing at the end of the armrest cushion. Mount the motor control on the side of the seat support with velcro (I).

Important: +/- poles must be correctly oriented. The motor control is damaged if it is installed inverted.