



# VELA Meywalk 2000



VELA Meywalk 2000 Large



VELA Meywalk 2000 Medium



VELA Meywalk 2000 Small

VELA Meywalk 2000 is a series of gait trainers designed to help stimulate and develop children and adults' walking pattern. It provides optimal support through all the phases of walking, encouraging mobility and independence. The ability to move freely and unaided is vital for the development of the social and interpersonal skills and for the user's health and comfort.

VELA Meywalk 2000 is available in 3 sizes and fit children from approximate six years and up to adulthood. The gait trainer has spring suspension allowing the seat and trunk support to follow the user's movements encouraging a better and more natural gait while providing optimal support through phases of walking.

The seat and trunk support of the gait trainer is easily adjusted to fit the user's size and level of upper body control. The seat unit can be lowered by using a swing bar, making it easier for the user to get safely in and out of the gait trainer and putting less strain on the caregivers.

Walking in a VELA Meywalk 2000, users are not dependent on their forearms for weightbearing which encourages social participation. Additional support and stability can be found in the integrated hand bar which can be adjusted to the user.

Standing in an upright position daily has a profound influence on the user's muscle strength and bone density. Developing the ability to stand and walk with a VELA Meywalk 2000 vastly improves blood circulation, digestion, and respiration which is often a necessity for managing other daily activities.

## BENEFITS

- :: User can be standing/walking despite the lack of or reduced standing function
- :: The spring suspension allows the seat and trunk support to follow the user's movements, promoting a natural walking pattern.
- :: Safe and easy to access with the unique lifting and lowering system of the seat.
- :: It is simple and easy to adjust.
- :: Grows with the user thanks to the variable height adjustment of seat and body ring

## STANDARD CONFIGURATION



VELA Meywalk 2000 Small

VELA Meywalk 2000 Medium

VELA Meywalk 2000 Large

MEASUREMENTS AND WEIGHT

Length	85 cm	95 cm	95 cm
Width, outer measures	65 cm	65 cm	65 cm
Width, inner measures	50 cm	50 cm	50 cm
Seat width:			
Front	9 cm	10 cm	10 cm
Middle	6 cm	7 cm	7 cm
Back	14 cm	15 cm	15 cm
Seat depth	27 cm	32 cm	32 cm
Seat height/stride length*	48-66 cm	62-86 cm	76-100 cm
Height of body ring*	87-103 cm	99-126 cm	116-143 cm
Size of body ring (select size based on waist measurement of the user)	70, 90, 105, 120 og 130 cm	70, 90, 105, 120 og 130 cm	70, 90, 105, 120 og 130 cm
Weight (without accessories)	22 kg	25 kg	26 kg
Maximum weight load	100 kg	100 kg	100 kg
Weight interval (colour on springs) (Select based on the body weight of the user)	Min 20 kg (white) Min. 45 kg. (yellow) Min. 65 kg (blue) Min 80 kg. (red)	Min 20 kg (white) Min. 45 kg. (yellow) Min. 65 kg (blue) Min 80 kg. (red)	Min 20 kg (white) Min. 45 kg. (yellow) Min. 65 kg (blue) Min 80 kg. (red)
Colour of the gait trainer	Green	Blue	Blue
Frame	Steel	Steel	Steel
Item number	MW104500	MW104175	MW104173

\* Height is measured with half compression of springs.

All VELA Meywalk 2000 gait trainers are CE marked and comply with current standards.

The gait trainer is suitable for use on horizontal, stable, and plain surfaces under controlled conditions.

The standard gait trainer is supplied with upholstered parts in artificial leather and polyurethane foam which are durable and easy to clean and disinfect.

### ACCESSORIES AND CUSTOMISED SOLUTIONS

All VELA Meywalk 2000 gait trainers are versatile and customisable to the user's specific needs thanks to a wide range of accessories, for example, a hip and trunk support to optimise the body posture. You can find a complete accessory list with item numbers, pictures, and descriptions at [vela.eu](http://vela.eu).

Distributor Details:

### WORK ENVIRONMENT

The VELA Meywalk 2000 gait trainer is versatile and adjusts easily without tools, especially beneficial in health-care environments where several users may share one gait trainer.

The seat can be lowered using the integrated swing bar so users can relocate to the gait trainer with minimal assistance. It is also possible to lift the user into the Meywalk 2000 gait trainer.

### INTENDED USE

Walker for children and adults, to help them stand up and move around despite having reduced or no walking function