
Vela Miniwalk



VELA Miniwalk is a gait trainer for children with mobility disabilities. It is available in two sizes for children at ages 1-7 years.

VELA Miniwalk is designed to encourage movement promoting the development of the child's gait pattern. It provides optimal support through all the phases of walking enabling the child to move and play independently.

The earlier the child is stimulated to walk and explore, the greater the chance of achieving a gait pattern.

The seat and trunk support of the VELA Miniwalk enable the child to move around safely and independently while building postural strength and stability. In addition, VELA Miniwalk is spring suspended, which allows the gait trainer to follow the child's movements stimulating the muscles helping the child develop a gait pattern.

VELA Miniwalk has a trunk support and seat that are both continuously adjustable with regards to height allowing the gait trainer to 'grow' with the child. These settings are easily carried out without the use of tools.

The seat unit can be lowered by using a swing bar, making it easier for the user to get safely in and out of the gait trainer and putting less strain on the caregivers.

The VELA Miniwalk is easy to mobilize, so the child can play and move around effortlessly without using the hands for weight bearing. If additional support is needed, the integrated and adjustable hand bar is positioned to help the child balance and gain control of the body.

BENEFITS

- :: The child can be standing/walking despite the lack of or reduced standing function.
- :: VELA Miniwalk moves up and down with each step stimulating the child's natural gait pattern.
- : Independent height adjustment of the trunk support and seat allows this gait trainer to grow with the child.
- :: Safe and easy to access using the integrated swing bar.
- :: Simple and easy to adjust.

STANDARD CONFIGURATION

	VELA Miniwalk Size 1	VELA Miniwalk Size 2
Length	74 cm	74 cm
Width, outer measures	59 cm	59 cm
Width, inner measures	42 cm	42 cm
Seat width: Front Middle Back Seat depth	7 cm 5 cm 11 cm 25 cm	7 cm 5 cm 11 cm 25 cm
Seat height/stride length*	27-48 cm	27-63 cm
Height of body ring*	60-72 cm	76-92 cm
Size of body ring (select size based on waist measurement of the user)	70 og 90 cm	70 og 90 cm
Weight (without accessories)	14 kg	14 kg
Maximum weight load	40 kg	40 kg
Weight interval (colour on springs) (Select based on the body weight of the user)	Min 10 kg (green) Min. 15 kg. (silver) Min. 20 kg (black)	Min 10 kg (green) Min. 15 kg. (silver) Min. 20 kg (black)
Colour of the gait trainer	Turquoise	Turquoise
Frame	Steel	Steel
Item number	MW106061	MW106062
* Height is massured with half compression of corings		

* Height is measured with half compression of springs.

All VELA Miniwalk gait trainers are CE marked and com- WORK ENVIRONMENT ply with current standards.

The gait trainer is suitable for use on horizontal, stable, and plain surfaces under controlled conditions.

The standard gait trainer is supplied with upholstered parts in artificial leather and polyurethane foam which are durable and easy to clean and disinfect.

ACCESSORIES AND CUSTOMISED SOLUTIONS

All VELA Miniwalk gait trainers are versatile and customisable to the child's specific needs thanks to a wide range of accessories, for example, a hip and trunk support to optimise the body posture. You can find a complete accessory list with item numbers, pictures, and descriptions at vela.eu

The VELA Miniwalk gait trainer is versatile and adjusts easily without tools, especially beneficial in healthcare environments where several children may share one gait trainer.

The seat can be lowered using the integrated swing bar so users can relocate to the gait trainer with minimal assistance. It is also possible to lift the user into the Miniwalk gait trainer.

INTENDED USE

Walker for children and adults, to help them stand up and move around despite having reduced or no walking function

Distributor Details: